


Insights for post-rejection care with Prospera™ dd-cfDNA surveillance¹



Key findings¹:

- **Post-rejection Prospera™ dd-cfDNA trends are prognostic of 1-year outcomes:** Longitudinal dd-cfDNA patterns following treatment (e.g. LOW, DROP vs MID, HIGH) were strongly associated with 1-year clinical outcomes, including renal function recovery, recurrent rejection, DSA formation and graft loss.¹
 - **97.6% (40/41)** of cases with sustained elevated dd-cfDNA levels (HIGH or MID) experienced negative outcomes.¹
 - Overall, patients with lower-risk dd-cfDNA trends were 60x more likely to have a positive outcome than those with higher risk trends.¹
- **Trends of dd-cfDNA post-rejection may serve as a molecular indicator of treatment effectiveness:** Prospera™ may help clinicians identify “favorable prognosis” versus “unfavorable prognosis” to immunosuppression before changes in serum creatinine levels.¹ In this study, creatinine trends had no statistical correlation with outcomes.²



“Following a rejection episode, understanding the long-term prognosis and risk is critical in making ongoing treatment decisions for the patient. These findings make it clear that dd-cfDNA is an important non-invasive tool to provide more accurate risk stratification post-rejection that will enable more personalized patient management and improve long-term graft survival.”

SUPHAMAI BUNNAPRADIST, M.D., PRINCIPAL INVESTIGATOR OF THE PEDAL STUDY²

Start the conversation about integrating Prospera™ into your surveillance protocol. **Contact Natera today.**



Read the publication

* Data based on this initial publication may be subject to change based on the final results of the PEDAL Study.

References

1. Bunnapradist, et al. Associations between donor-derived cell-free DNA dynamics and clinical outcomes after kidney allograft rejection: A prospective, multicenter study. Am. J. Transplant, 2025.
2. Natera, Inc. (2025, July 24). PEDAL Study Successful, shows monitoring with Prospera™ Kidney provides accurate prognosis of long-term clinical outcomes following rejection; now published in AJT. Business

Wire. Retrieved August 14, 2025, from <https://investor.natera.com/news/news-details/2025/PEDAL-Study-Successful-Shows-Monitoring-with-Prospera-Kidney-Provides-Accurate-Prognosis-of-Long-Term-Clinical-Outcomes-Following-Rejection-Now-Published-in-AJT/default.aspx>

13011 McCallen Pass, Building A Suite 100 | Austin, TX 78753 | www.natera.com

Prospera™ has been developed and its performance characteristics determined by the CLIA-certified laboratory performing the test. The test has not been cleared or approved by the US Food and Drug Administration (FDA). CAP accredited, ISO 13485 certified, and CLIA certified. © 2025 Natera, Inc. All Rights Reserved. PRO_CS_BR_PEDAL-Booklet_20250819_NAT-8022758



Prospera™ in practice: Prognostic power & insights from the PEDAL study

Real-World Evidence in Post-Rejection Monitoring Using Donor-derived Cell-Free DNA (dd-cfDNA)

The PEDAL Study is the **first prospective, multicenter clinical study** to evaluate how dd-cfDNA dynamics after biopsy-proven rejection can predict long-term kidney transplant outcomes.¹ Across 28 participating centers, the study demonstrated that Prospera™ dd-cfDNA trends over the 8-week post-rejection window were significantly associated with clinical outcomes at 12 months.¹

Four distinct dd-cfDNA patterns—categorized as **LOW**, **DROP**, **MID**, and **HIGH**—helped differentiate those with “favorable prognosis” from “unfavorable prognosis” to rejection treatment, providing a non-invasive molecular lens into post-treatment graft health.¹

This booklet highlights four case studies, each aligned with a different dd-cfDNA trajectory observed in the PEDAL Study¹:

Case 1: LOW

A patient with persistently low dd-cfDNA levels and full renal function recovery

Case 2: DROP

A case with initially high but sharply declining dd-cfDNA, associated with a favorable prognosis

Case 3: MID

A patient whose dd-cfDNA levels decreased then rebounded, correlating with poor outcome

Case 4: HIGH


A patient with persistently elevated dd-cfDNA and multiple adverse events

Together, these real-world cases underscore the clinical value of Prospera™ in the post-rejection setting, helping clinicians personalize management decisions, adjust surveillance intensity, and act proactively in response to molecular trends.

When rejection resolves: LOW and DROP trend profiles¹

CASE 1 (LOW Group)¹:

A patient with persistently low dd-cfDNA levels and full renal function recovery



- **Rejection Type:** T-cell mediated rejection (TCMR) 1A
- **Serum Creatinine:** Decrease from 2.48 mg/dL to 1.96 mg/dL




Outcome: No adverse events at 1 year



Key takeaway: Consistently low Prospera™ dd-cfDNA trends post-rejection may be associated with full recovery.¹

CASE 2 (DROP Group)¹:

A case with initially high but sharply declining dd-cfDNA, associated with a favorable prognosis



- **Rejection Type:** TCMR 1A
- **Serum Creatinine:** Decrease from 1.14 mg/dL to 1.06 mg/dL



Outcome: No adverse events at 1 year




Key takeaway: High initial Prospera™ dd-cfDNA levels that drop rapidly may reflect sufficient treatment response.¹

When rejection persists: MID and HIGH trend profiles¹

CASE 3 (MID Group)¹:

A patient whose dd-cfDNA levels decreased then rebounded, correlating with poor outcome



- **Rejection Type:** TCMR 1B progressing to chronic active TCMR 2




Outcome: Developed antibody-mediated rejection (ABMR) 4 months later



Key takeaway: A rebound in Prospera™ dd-cfDNA levels may indicate unresolved injury, despite initial improvement.¹

CASE 4 (HIGH Group)¹:

A patient with persistently elevated dd-cfDNA and multiple adverse events



- **Rejection Type:** ABMR
- **Serum Creatinine:** Increase from 1.5 mg/dL to 5.7 mg/dL



Outcome: Graft loss following second rejection biopsy



Key takeaway: Persistently elevated Prospera™ dd-cfDNA may be a strong signal of treatment resistance.¹